



COMPUTER DRIVEN SKIN ANALYSIS EQUIPMENT

These devices require a computer in order to run. they give the possibility of saving results and in some case providing a product or treatment recommendation.

[Dermacheck 3000 Software](#)

DERMACHECK 3000 SOFTWARE

The Dermacheck 3000 software offers aesthetic and dermatology clinics additional measurement modules for use with a selection of skin analysis equipment.

[READ MORE](#)

[Multi Skin Test Center MC 1000](#)

MULTI SKIN TEST CENTER MC 1000

The sleek all white housing and white probes of the Multi Skin Test Center MC 1000 will add a touch of class to your skin assessment with your client.

[READ MORE](#)

VISIOFACE® 1000 D

The Visioface® is a full face photographic method for facial skin analysis, optimal product recommendation and treatment documentation. It consists of a booth in which the face is placed for a photographic image.

READ MORE

VISIOFACE® LITE

The Visioface® Lite is a full face photographic method for facial skin analysis, optimal product recommendation and treatment documentation. It consists of a small lightweight booth made out of a special light transmitting plexiglass in which the face is placed for a photographic image.

READ MORE

VISIOPOR® PP 34N

Acne is a common disorder of the pilosebaceous follicles with the multifactorial pathogenesis. It typically begins in adolescence when androgen hormones stimulate the production of sebum and proliferation of follicular epidermis.

READ MORE