



SKIN MEASUREMENT PARAMETERS

A wide variety of skin parameters can be measured.

Cutometer® Dual MPA 580

CUTOMETER® DUAL MPA 580

The mechanical properties are important physiologically, socially and emotionally. Physiologically they provide the strength and flexibility of the skin.

READ MORE

Skin-Thermometer ST 500

SKIN-THERMOMETER ST 500

Skin temperature measurement is a valuable tool in many fields of application. It can be used in efficacy testing, claim support, objective clinical analysis, occupational medicine, medical consultancy and many more fields.

READ MORE