



MULTIPROBE ADAPTER SYSTEM

Courage & Khazaka

BENEFITS

- ⦿ Expand your measurement capabilities easily
- ⦿ Cost effective
- ⦿ Flexible
- ⦿ New developments compatible with system
- ⦿ User friendly software for speedy transfer of data for analysis
- ⦿ Calibration data stored in the probe
- ⦿ Large range of measurement parameters available
- ⦿ Check calibration functions to check accuracy
- ⦿ Basic systems to suit different needs

The MultiProbe Adapter System® provides the user with increased flexibility and cost effectiveness for multiple measurements. Each system has the capacity to connect several probes of either the same or different measurements. All of the systems require connection to a computer and come with the software that will run the device. The software for the MPA system automatically recognizes the probes that are connected to the system making it quick and easy to start making measurements. In addition all future Courage & Khazaka measurement probes will be compatible with this system ensuring that you can expand your measurement capabilities in the most cost effective way.

MEASUREMENTS AVAILABLE

- ⦿ Hydration ~ Corneometer® CM 825
- ⦿ Sebum ~ Sebumeter® SM 815
- ⦿ Pigmentation ~ Mexameter® MX 18
- ⦿ Skin Colour ~ Skin-Colorimeter CL 400
- ⦿ Skin Gloss ~ Skin-Glossometer GL 200

- Ambient Condition Sensor ~ RHT 100
- Elasticity ~ Cutometer® dual MPA 580
- Skin pH ~ Skin-pH-Meter® PH 905
- Skin Temperature ~ Skin-Thermometer ST 500
- Direction of Fibres ~ Reviscometer® RVM 600
- Skin Friction ~ Frictiometer FR 700
- Transepidermal Water Loss (TEWL) ~ Tewameter® TM 300 & Tewameter® Triple TM 330T

WHICH BASIC DEVICES ARE AVAILABLE?

MULTIPROBE ADAPTER MPA 6

- Inbuilt Sebumeter® socket
- Connection of five probes (same or different)
- Operated using MPA software
- Can accept all probes except Cutometer® probes

MULTIPROBE ADAPTER MPA 10

- Inbuilt Sebumeter® socket
- Connection of ten probes (same or different)
- Operated using MPA software
- Can accept all probes except Cutometer® probes

MULTIPROBE ADAPTER MPA 2

- Small USB powered
- Connection of two probes
- Operated using MPA software
- Can accept all probes except Cutometer® probes

CUTOMETER® DUAL MPA 580

- Inbuilt Sebumeter® socket
- Connect two Cutometer® probes (different diameters)

- Connection of four additional probes
- Operated with Cutometer® dual software and MPA software

MULTI DISPLAY DEVICE MDD 4

- Stand alone with one probe and RHT 100
- Connection of two additional probes
- Large Colour display for showing results
- Can also be operated using MPA software